

FABRICATION OF A NEW PROSTHESIS WITH PRECI-BAR

Article numbers:	1100/H/MR/OR	1100/H/MR/PA	1100/H/MR/IN
	1102/H/MR/OR	1102/H/MR/PA	1102/H/MR/IN
	1103/B/OR	1103/B/PA	1012/B
	1103/BL/OR	1103/BL/PA	
	1101/B/OR	1101/B/PA	1013/B
	1101/BL/OR	1101/BL/PA	
	1100/B/OR	1100/B/PA	1011/B
	1100/BL/OR	1100/BL/PA	

1. Take a full arch impression.
2. Pour a stone model. The entire bar area must be made of epoxy material or self-curing acrylic resin.
3. Reduce the corresponding rider to the required length (according to the bar) (**1100/H/MR/OR, 1100/H/MR/PA, 1100/H/MR/IN, 1102/H/MR/OR, 1102/H/MR/PA, 1102/H/MR/IN**).
4. Place the rider on the model. When processing a resilient PRECI-BAR, the space maintainer must be used.
5. Prepare the wax setup for investing.
6. Separate the reproduced bar from the acrylic resin. Complete the prosthesis and remove the space maintainer if used.
7. Clean the rider and activate it if necessary (see INFO 108).