

## SECURING OF SPRING PINS

As a result of occlusal forces, spring pins may gradually loosen from base rings or retention parts. Therefore, the spring pins should be secured with **CEKA BOND (CB1)** to prevent spontaneous loosening.

1. Clean the threads of the spring pin and base ring or retention part prior to screwing in.
2. Apply a small amount of **CEKA BOND** to the threads of the spring pin.
3. Screw the spring pin with the **A 1** or **RE H 5** accessory fully down into the base ring or the retention part (see INFO 063).
4. Remove the excess material.
5. **CEKA BOND** is not a permanent fixation. With the **A 1** or **RE H 5** accessory, the spring pin may still be easily unscrewed (see INFO 063).
  - **CEKA BOND** is an anaerobic adhesive.
  - **Safety data sheet:** available on request from your local CEKA distributor.
  - **Shelf-life:** see date on package.
  - **Store at room temperature.**