

NEW PROSTHESIS WITH DAMAGED FEMALES (M2)

If a new prosthesis must be fabricated in case of abutment retainers with damaged females, the following procedure must be followed.

CHAIRSIDE PROCEDURE

1. Assemble the **RE H 14** impression tool with the **RE 0031/195** oversized spring pin.
2. Place this combination in the female in the mouth and use the **RE 0096** space maintainer.
3. The combination must be absolutely stable and secure in the female. If not, activate the spring pin (see INFO 059).
4. Block out the undercuts of the female with soft wax.
5. Take impression with a custom-made impression tray.

LABORATORY PROCEDURE

1. Assemble the **RE H 14** impression tool (including the **RE 0031** normal spring pin) with the **RE 0096** space maintainer and the **RE H 13** model accessory.
2. Place this combination in the impression.
3. Pour a model in stone which incorporates the model analogues.
4. Always use the **RE 0096** space maintainer for further processing.
5. Make a new secondary construction.
6. After fabrication of a new prosthesis, replace the **RE 0031** normal spring pins with the oversized spring pins used during impression-taking and secure them with **CEKA BOND** (see INFO 069).